



Kinetic Profiling by NMR

This note summarises a procedure for obtaining kinetic (time course) data with NMR on **Bruker** instruments using **TOPSPIN**. This provides a means of monitoring reactions as a function of time in the NMR tube and potentially for determining kinetic parameters from the time course data.

1. First tune, shim and run initial spectra of your starting materials **before** the reaction to make sure you know exactly where they appear in the spectrum with the same solvent that you use in the reaction under the same conditions, e.g. temperature, preferably on the same day!
2. Having optimised the spectra of your starting materials alone (tuning the probe and shimming), quickly run the first spectrum of your reaction after mixing noting the time you started (the command **dpa** will display the exact end time of each acquisition and these times can be used in generating time course profiles). If your reaction is not too fast, you may have time to do some quick shimming beforehand; re-tuning is not necessarily required.
3. Having run the first reaction spectrum, increment experiment number (**iexp**), and use the command **multi_zgvd** (*answer the following questions*).

Enter fixed delay (in seconds)

e.g. **600** (10 minutes for spectrometer to wait until next acquisition)

Number of Experiments

e.g **30**

The program should tell you how long your kinetic run should take.

4. Process your first spectrum using **efp** and phase it perfectly, manually. You need all subsequent spectra in the kinetic run to be phased in exactly the same way. Use the **multicmd** command to do this, which asks you the number of experiments required, and the command for processing. Execute with the **efp** command then repeat using the **abs** command (for baseline correction). The macro will use [**efp**] then [**abs**] for all the datasets you have specified; it will phase the specified spectra exactly the same as the first one you phased perfectly.
5. Having now obtained all your experiments, please look closely for your signals of interest, those that grow at the expense of those that decrease. You will then need to integrate them in exactly the same way e.g. integral limits and peak positions, for all *n* experiments.
6. **Integrate** a good example spectrum of the kinetic run, where you have **both** starting material and product and **save** regions to '*intrng*'.

Then write a new integral file with command **wmisc (write new)** with filename xxx, enter a new name.

To be on the safe side, check the integrals are what you want with **rmisc** (read xxx) on another spectrum within your list.

7. When satisfied with your integration of your starting material and product peaks, then go to the first spectrum of the reaction mixture and use the command **multi_integ3** (answer list of questions);

Use EXPNOs(0) or PROCNOs((1)

0

Enter first experiment No

e.g. **2**

Enter no of experiments

e.g. **30**

Enter name of intrng file**xxx**


The program calculates all the integrals in files 1.....n and the results are conveniently presented as a txt file in the experiment you performed the **multi_integ3** command in, with the processed spectrum e.g. D:/NMR/data/organic.nmr/file/pdata/1/file_int.txt.

An example of how this text file should appear is shown below;

The screenshot shows a Notepad window titled "kinetics 15 mol% ligand_int - Notepad". The window contains a table of numerical data with 21 columns and 21 rows. The data is organized as follows:

Row	Col 1	Col 2	Col 3	Col 4	Col 5	Col 6	Col 7	Col 8	Col 9	Col 10	Col 11	Col 12	Col 13	Col 14	Col 15	Col 16	Col 17	Col 18	Col 19	Col 20	Col 21
1	0.000000	332703000.000000	1806560.000000	17015600.000000	48085200.000000	263676000.000000															
2	0.000000	333041000.000000	2207930.000000	17206000.000000	48128600.000000	263505000.000000															
3	0.000000	334120000.000000	2656190.000000	17080800.000000	48161800.000000	262826000.000000															
4	0.000000	335279000.000000	3131670.000000	17223900.000000	48093000.000000	262569000.000000															
5	0.000000	336181000.000000	3647280.000000	17112000.000000	47965500.000000	261604000.000000															
6	0.000000	337178000.000000	4305880.000000	17158500.000000	48020600.000000	261266000.000000															
7	0.000000	338341000.000000	4902940.000000	17151900.000000	47852200.000000	260211000.000000															
8	0.000000	339662000.000000	5475170.000000	17215900.000000	47781700.000000	259940000.000000															
9	0.000000	340479000.000000	5998760.000000	17196700.000000	47802100.000000	259056000.000000															
10	0.000000	341304000.000000	6661750.000000	17225400.000000	47729600.000000	258699000.000000															
11	0.000000	342899000.000000	7313360.000000	17227300.000000	47811000.000000	257917000.000000															
12	0.000000	343503000.000000	7958970.000000	17232500.000000	47586400.000000	257550000.000000															
13	0.000000	345011000.000000	8567460.000000	17207900.000000	47751800.000000	256719000.000000															
14	0.000000	345761000.000000	9278720.000000	17323100.000000	47664100.000000	256249000.000000															
15	0.000000	346809000.000000	9921190.000000	17262200.000000	47536600.000000	255075000.000000															
16	0.000000	347439000.000000	10822600.000000	17381800.000000	47471300.000000	254891000.000000															
17	0.000000	349043000.000000	11526600.000000	17311500.000000	47409900.000000	253482000.000000															
18	0.000000	350172000.000000	12286600.000000	17348100.000000	47499700.000000	253192000.000000															
19	0.000000	351138000.000000	12942300.000000	17269000.000000	47392900.000000	252026000.000000															
20	0.000000	351750000.000000	13646800.000000	17370100.000000	47415000.000000	251893000.000000															
21	0.000000	352524000.000000	14343100.000000	17329200.000000	47266100.000000	250391000.000000															
22	0.000000	352456000.000000	15240300.000000	17326400.000000	47251700.000000	250089000.000000															

This table of integral values can easily be imported into **Excel** as a Space and Tab delimited file for plotting and for calculations of rate constants, as shown below.



The in-house au program **readdate** may also be used to generate a text file containing the date and time points for a series of spectra of increasing experiment number (expno) and avoids the need to manually extract the time from each spectrum with the dpa command. The generated file (*kineticdates.txt*) will be stored in the expno of the first experiment. Note: this is not part of the standard Bruker installation and will have to be added to your TOPSPIN installation: copy the au program into the directory:

C:\Bruker\TOPSPIN\exp\stan\nmr\au\src\

(or your corresponding Bruker installation path) and, in your first experiment of the kinetic run, use the command **xau readdate** to compile and run the program.

Please see Tim Claridge or Nader Amin if you have any problems with these procedures.